

Product Datasheet

CD15 Antibody (Biotin) (orb1810979)

Description	CD15 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC, IHC-P
Immunogen	Mouse F9 Teratocarcinoma Stem Cells (X-irradiated)
Target	CD15
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Concentration	1mg/ml
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse Ig1, k
Clonality	Monoclonal
Uniprot ID	/
Expiration Date	12 months from date of receipt.