

Product Datasheet

Tenovin-6 (orb181161)

Description	Tenovin-6 is the water soluble analog of Tenovin-1 and acts as a potent SIRT1 (IC50=21 uM) and SIRT2...
Note	For research use only
Formula	$C_{25}H_{34}N_4O_2S$
MW	454.6
CAS Number	[1011557-82-6]
SMILES	<chem>CC(C)(C)C1=CC=C(C=C1)C(=O)NC(=S)NC2=CC=C(C=C2)NC(=O)CCCCN(C)C</chem>
Expiration Date	12 months from date of receipt.