

Product Datasheet

TGX221 (orb181167)

Description	TGX-221 is a potent, selective, and cell permeable inhibitor of PI3K p110 β .
Note	For research use only
Formula	C ₂₁ H ₂₄ N ₄ O ₂
MW	364.44
CAS Number	[663619-89-4]
SMILES	<chem>O=C1N2C(C(C(NC3=CC=CC=C3)C)=CC(C)=C2)=NC(N4CCOCC4)=C1</chem>
Expiration Date	12 months from date of receipt.