

Product Datasheet

AMH Antibody (HRP) (orb1811755)

Description	AMH Antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	HRP
Tested Applications	IHC-P, WB
Immunogen	Recombinant human AMH protein (25-450/560aa)
Target	AMH
Concentration	1mg/ml
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P03971
Expiration Date	12 months from date of receipt.