

Product Datasheet

AMH Antibody (Biotin) (orb1811756)

Description	AMH Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	IHC-P, WB
Immunogen	Recombinant human AMH protein (25-450/560aa)
Target	AMH
Concentration	1mg/ml
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P03971
Expiration Date	12 months from date of receipt.