

Product Datasheet

THY1 Antibody (PerCP/Cy5.5) (orb1811974)

Description	THY1 Antibody (PerCP/Cy5.5)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	Recombinant Human Thy-1 protein, C-His (HEK293) (20-141/161aa)
Target	THY1
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P04216
Expiration Date	12 months from date of receipt.