

Product Datasheet

THY1 Antibody (HRP) (orb1811991)

Description	THY1 Antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-P, WB
Immunogen	Recombinant Human Thy-1 protein, C-His (HEK293) (20-141/161aa)
Target	THY1
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P04216
Expiration Date	12 months from date of receipt.