

Product Datasheet

HSPA8 Antibody (APC/Cy5) (orb1812050)

Description	HSPA8 Antibody (APC/Cy5)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	APC/Cy5
Tested Applications	ICC
Immunogen	Recombinant human HSC70 Protein
Target	HSPA8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P11142
Expiration Date	12 months from date of receipt.