


Product Datasheet

Fibroblast growth factor 8 FGF8 Antibody (orb18123)

Description

Fibroblast growth factor 8 FGF8 Antibody



Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunohistochemical analysis of paraffin...

Immunogen


A synthetic peptide corresponding to a sequence at the C-terminus of human FGF8 (163-185aa FMKRLPRGHHTTEQSLRFEFLNY), identical to the related rat and mouse sequences.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.



Lane 1 : MCF7 Whole Cell Lysate
Lane 2 : HeLa Whole Cell Lysate
Lane 3 : SMMC Whole Cell Lysate
Lane 4 : HT1080 Whole Cell Lysate
Lane 5 : colo320 Whole Cell Lysate

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Western Blot analysis of cell lysate usi...

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

26525 MW

Uniprot ID

P55075

Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Rat, Mouse, By Heat Western blot, 0.1-0.5µg/ml, Human, Rat, Mouse

Expiration Date

12 months from date of receipt.