

Product Datasheet

TMPRSS2 Antibody (Biotin) (orb1812499)

Description	TMPRSS2 Antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant human TMPRSS2 protein (108-492/492aa)
Target	TMPRSS2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15393
Expiration Date	12 months from date of receipt.