

Product Datasheet

PIK3C3 Antibody (APC/Cy7) (orb1812540)

Description	PIK3C3 Antibody (APC/Cy7)
Species/Host	Rabbit
Reactivity	Human
Conjugation	APC/Cy7
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PIK3C3 around the phosphorylation site of Ser164 RT(p-S)ST
Target	PIK3C3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NEB9
Expiration Date	12 months from date of receipt.