

Product Datasheet

EEA1 Antibody (PerCP/Cy5.5) (orb1812660)

Description	EEA1 Antibody (PerCP/Cy5.5)
Species/Host	Rabbit
Reactivity	Hamster, Human
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human EEA1 (1251-1350/1411aa)
Target	EEA1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15075
Dilution Range	Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.