

Product Datasheet

GFP Antibody (AP) (orb1812690)

Description	GFP Antibody (AP)
Species/Host	Mouse
Conjugation	AP
Tested Applications	WB
Immunogen	Recombinant Green Fluorescence Protein
Target	GFP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P42212
Expiration Date	12 months from date of receipt.