

Product Datasheet

ENG Antibody (PE) (orb1812758)

Description	ENG Antibody (PE)
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	Recombinant human CD105 (26-586/658 (C-6x His-Tag)aa)
Target	ENG
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P17813
Expiration Date	12 months from date of receipt.