

Product Datasheet

SCNN1B Antibody (PerCP/Cy7) (orb1812940)

Description	SCNN1B Antibody (PerCP/Cy7)
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse
Conjugation	PerCP/Cy7
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human SCNN1B around the phosphorylation site of Ser633 IE(p-S)D
Target	SCNN1B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51168
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.