

Product Datasheet

ENPP1 Antibody (Biotin) (orb1812990)

Description	ENPP1 Antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from mouse ENPP1 (851-900/925aa)
Target	ENPP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P22413
Expiration Date	12 months from date of receipt.