

Product Datasheet

RBP4 Antibody (Biotin) (orb1813338)

Description	RBP4 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant human RBP4 protein (19-201/201aa)
Target	RBP4
Preservatives	0.01M PBS(pH7.4)
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P02753
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.