

Product Datasheet

TP63 Antibody (APC/Cy5.5) (orb1813421)

Description	TP63 Antibody (APC/Cy5.5)
Species/Host	Mouse
Reactivity	Human
Conjugation	APC/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human P63
Target	TP63
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9H3D4
Expiration Date	12 months from date of receipt.