

Product Datasheet

HLA-C Antibody (AP) (orb1813552)

Description	HLA-C Antibody (AP)
Species/Host	Mouse
Reactivity	Human
Conjugation	AP
Tested Applications	IHC-P
Immunogen	Recombinant human HLA-C protein (25-305/366aa)
Target	HLA-C
Storage	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P10321
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.