

Product Datasheet

ML347 (orb181404)

Description	ML347 is a selective BMP receptor inhibitor with IC50 of 32 nM for ALK2, > 300-fold selectivity over...
Note	For research use only
Formula	C ₂₂ H ₁₆ N ₄ O
MW	352.39
CAS Number	[1062368-49-3]
SMILES	<chem>N1C2C(=C(C3=C4N(N=C3)C=C(C3=CC=C(OC)C=C3)C=N4)C=CC=2)C=CC=1</chem>
Expiration Date	12 months from date of receipt.