


## Product Datasheet

### DBeQ (orb181411)

|                        |                                                                                                                                                                 |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | DBeQ is a selective, potent, reversible, and ATP-competitive p97 inhibitor with IC50 of 1.5 µM.                                                                 |
| <b>Note</b>            | For research use only                                                                                                                                           |
| <b>Formula</b>         | C <sub>22</sub> H <sub>20</sub> N <sub>4</sub>                                                                                                                  |
| <b>MW</b>              | 340.421                                                                                                                                                         |
| <b>CAS Number</b>      | [177355-84-9]                                                                                                                                                   |
| <b>SMILES</b>          | <chem>N([H])(C([H])([H])C1C([H])=C([H])C([H])=C([H])C=1[H])C1C2=C([H])C([H])=C([H])C([H])=C2N=C(N=1)N([H])C([H])([H])C1C([H])=C([H])C([H])=C([H])C=1[H])</chem> |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                                                 |



Chemical structure of DBeQ