


Product Datasheet

CLEC12A antibody (orb181425)

Description	Goat polyclonal antibody to CLEC12A.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CLEC12A / MICL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize all reported isoforms
Protein Sequence	NPVQLGSTYFRE
Clonality	Polyclonal
MW	30.8; 27.3; 32.0
Entrez	160364
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of U937 lysate using CLEC12A antibody.