


Product Datasheet

NFATC3 antibody (orb181433)

Description	Goat polyclonal antibody to NFATC3.
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NFATC3 / NFAT4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 150kDa band observed in nuclear lysates of cell line Jurkat (calculated MW of 116kDa according to NP_775188.1). This molecular weight is routinely observed by other sources. Recommended concentration: 0.03-0.1µg/ml.
Protein Sequence	HQELDAHEDDLQIN
Clonality	Polyclonal
MW	116; 115; 115; 113
Entrez	4775
Dilution Range	ELISA: 1:16000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat nuclear ...