

Product Datasheet

RRM2B Antibody (AR700) (orb1814368)

Description	RRM2B Antibody (AR700)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	AR700
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human p53R2/RRM2B (131-230/351aa)
Target	RRM2B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7LG56
Expiration Date	12 months from date of receipt.