

Product Datasheet

ABCB1 Antibody (Biotin) (orb1814525)

Description	ABCB1 Antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	Recombinant Human P-glycoprotein 1 protein
Target	ABCB1
Preservatives	1*TBS (pH7.4), 0.05% rAlbumin, 40% Glycerol. Preservative: 0.05% Sodium Azide.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P08183
Expiration Date	12 months from date of receipt.