

Product Datasheet

MUC1 Antibody (HRP) (orb1814759)

Description	MUC1 Antibody (HRP)
Species/Host	Mouse
Reactivity	Human
Conjugation	HRP
Tested Applications	IHC-P
Immunogen	Recombinant Human Mucin-1 Protein (921-1158/1255aa)
Target	MUC1
Preservatives	0.01M PBS(pH7.4)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P15941
Expiration Date	12 months from date of receipt.