

Product Datasheet

KRT14 Antibody (HRP) (orb1814787)

Description	KRT14 Antibody (HRP)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CK14
Target	KRT14
Preservatives	PBS, Glycerol, rAlbumin
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P08779
Expiration Date	12 months from date of receipt.