

Product Datasheet

3326-32-7 Antibody (Biotin) (orb1814995)

Description	3326-32-7 Antibody (Biotin)
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	KLH conjugated FITC
Target	3326-32-7
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Storage: Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile distilled water or diluent supplied, the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	3326-32-7
Expiration Date	12 months from date of receipt.