

Product Datasheet

LAPTM5 Antibody (APC/Cy5) (orb1815071)

Description	LAPTM5 Antibody (APC/Cy5)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	APC/Cy5
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human LAPTM5 (101-200/262aa)
Target	LAPTM5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13571
Expiration Date	12 months from date of receipt.