

Product Datasheet

MKI67 Antibody (Biotin) (orb1815095)

Description	MKI67 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	Recombinant human Ki-67
Target	MKI67
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P46013
Expiration Date	12 months from date of receipt.