

Product Datasheet

VCAM1 Antibody (AP) [Out of stock] (orb1815115)

Description	VCAM1 Antibody (AP) [Out of stock]
Species/Host	Rabbit
Conjugation	AP
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human VCAM-1 (651-739/739aa)
Target	VCAM1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P19320
Expiration Date	12 months from date of receipt.