

Product Datasheet

MAPK7 Antibody (PerCP/Cy5.5) (orb1815295)

Description	MAPK7 Antibody (PerCP/Cy5.5)
Species/Host	Rabbit
Reactivity	Human
Conjugation	PerCP/Cy5.5
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ERK5 around the phosphorylation site of Ser496 GP(p-S)AP
Target	MAPK7
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13164
Expiration Date	12 months from date of receipt.