

Product Datasheet

TFPI2 Antibody (orb18153)

Description

TFPI2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence in the middle region of human TFPI2 (170-186aa NVTRYFNPRTYRTCDAF).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

26934 MW

Uniprot ID

P48307

Sensitivity



> 5000 cells

Dilution Range


Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By HeatWestern blot, 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Lane 1 : MM453 Whole Cell Lysate
Lane 2 : MM231 Whole Cell Lysate
Lane 3 : HeLa Whole Cell Lysate
Lane 4 : HT1080 Whole Cell Lysate
Lane 5 : Jurkat Whole Cell Lysate

Western Blot analysis of cell tissue lys...