

Product Datasheet

RBFOX3 Antibody (Biotin) (orb1815350)

Description	RBFOX3 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Mouse, Rat
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	Recombined human NeuN (1-100/312aa)
Target	RBFOX3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	A6NFN3
Expiration Date	12 months from date of receipt.