

Product Datasheet

CD274 Antibody (APC) (orb1815438)

Description	CD274 Antibody (APC)
Species/Host	Mouse
Conjugation	APC
Tested Applications	FC
Immunogen	Recombinant human CD274 protein
Target	CD274
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300.
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q9NZQ7
Expiration Date	12 months from date of receipt.