

Product Datasheet

NRXN2 Antibody (HRP) (orb1815538)

Description	NRXN2 Antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Neurexin 2 (231-330/1712aa)
Target	NRXN2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9P2S2
Expiration Date	12 months from date of receipt.