

Product Datasheet

PGC Antibody (Biotin) (orb1815550)

Description	PGC Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant human PGC protein (17-388/388aa)
Target	PGC
Preservatives	PBS (ph=7.4)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P20142
Expiration Date	12 months from date of receipt.