

Product Datasheet

LRP1 Antibody (BF350) (orb1815568)

Description	LRP1 Antibody (BF350)
Species/Host	Rabbit
Reactivity	Human
Conjugation	BF350
Tested Applications	ICC
Immunogen	KLH conjugated synthetic peptide derived from human LRP1 (2501-2600/4544aa)
Target	LRP1
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q07954
Expiration Date	12 months from date of receipt.