

Product Datasheet

ITGAL Antibody (APC/Cy5.5) (orb1815798)

Description	ITGAL Antibody (APC/Cy5.5)
Species/Host	Mouse
Reactivity	Human
Conjugation	APC/Cy5.5
Tested Applications	FC, ICC
Target	ITGAL
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2b, k
Clonality	Monoclonal
Uniprot ID	P20701
Dilution Range	Flow-Cyt=10ul/Test
Expiration Date	12 months from date of receipt.