

Product Datasheet

RPS6KB1 Antibody (AR700) (orb1815923)

Description	RPS6KB1 Antibody (AR700)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	AR700
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RPS6KB1 around the phosphorylation site of Ser434 IR(p-S)PR
Target	RPS6KB1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P23443
Expiration Date	12 months from date of receipt.