

Product Datasheet

IL1B Antibody (Biotin) (orb1815983)

Description	IL1B Antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Biotin
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse IL-1beta (1-100/269aa)
Target	IL1B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P01584
Expiration Date	12 months from date of receipt.