

Product Datasheet

STXBP1 Antibody (HRP) (orb1815989)

Description	STXBP1 Antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Munc18-1 around the phosphorylation site of Thr574 (501-594/594aa)
Target	STXBP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P61764
Expiration Date	12 months from date of receipt.