

Product Datasheet

LPAM-1 Antibody (PE) (orb1816061)

Description	LPAM-1 Antibody (PE)
Species/Host	Rat
Conjugation	PE
Tested Applications	FC
Target	LPAM-1
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	/
Expiration Date	12 months from date of receipt.