

Product Datasheet

Ly49I Antibody (FITC) (orb1816063)

Description	Ly49I Antibody (FITC)
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Target	Ly49I
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Clonality	Monoclonal
MW	4 kDa
Uniprot ID	O88358, O88357
Expiration Date	12 months from date of receipt.