

Product Datasheet

Thy1.2 Antibody (Pacific Blue) (orb1816066)

Description	Thy1.2 Antibody (Pacific Blue)
Species/Host	Rat
Conjugation	Pacific Blue
Tested Applications	FC
Target	Thy1.2
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	/
Expiration Date	12 months from date of receipt.