

## Product Datasheet

### CD72 Antibody (PE) (orb1816068)

|                            |                                                     |
|----------------------------|-----------------------------------------------------|
| <b>Description</b>         | CD72 Antibody (PE)                                  |
| <b>Species/Host</b>        | Mouse                                               |
| <b>Conjugation</b>         | PE                                                  |
| <b>Tested Applications</b> | FC                                                  |
| <b>Target</b>              | CD72                                                |
| <b>Preservatives</b>       | 0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300           |
| <b>Storage</b>             | Store at 2-8°C. Protect from light. Avoid freezing. |
| <b>Note</b>                | For research use only                               |
| <b>Clonality</b>           | Monoclonal                                          |
| <b>MW</b>                  | 40 kDa                                              |
| <b>Uniprot ID</b>          | <b>P21854</b>                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                     |