

Product Datasheet

EPCAM Antibody (Biotin) (orb1816080)

Description	EPCAM Antibody (Biotin)
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	FC
Target	EPCAM
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Clonality	Monoclonal
MW	40 kDa
Uniprot ID	Q99JW5
Expiration Date	12 months from date of receipt.