

Product Datasheet

CD14 Antibody (Biotin) (orb1816084)

Description	CD14 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC
Target	CD14
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Clonality	Monoclonal
MW	53-55 kDa
Uniprot ID	P08571
Expiration Date	12 months from date of receipt.