

Product Datasheet

ITGB1 Antibody (Biotin) (orb1816087)

Description	ITGB1 Antibody (Biotin)
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	FC
Target	ITGB1
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Clonality	Monoclonal
MW	130 kDa
Uniprot ID	P05556
Expiration Date	12 months from date of receipt.