

Product Datasheet

SARS-CoV-2 Antibody (Biotin) (orb1816137)

Description	SARS-CoV-2 Antibody (Biotin)
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein full length
Target	SARS-CoV-2
Preservatives	1*PBS, 0.03% Proclin300.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.